

September 10, 2014

OG Tube/Bougie VS. Suction Calibration System During Laparoscopic Sleeve Gastrectomy

**Rose Huang, MSBME
Michel Gagner, MD, FRCSC, FACS, FASMBS**

Hospital Du Sacre Coeur, Montreal, QC, Canada

*August 26-30, 2014
Montreal, Canada*



Boehringer Laboratories, LLC. • 300 Thoms Dr., • Phoenixville, PA 19460
1-800-642-4945 • www.visigi.com

Clinical Education Presented By:

BOEHRINGER®

OG Tube/Bougie vs. Suction Calibration System During Laparoscopic Sleeve Gastrectomy

Rose Huang, MSBME
Michel Gagner, MD, FRCSC, FACS, FASMBS
Hospital Du Sacre Coeur, Montreal, QC, Canada

Introduction


Expert consensus recommends a bougie should be in place when a sleeve gastrectomy is fashioned. However, this requires three gastric insertions: OG tube for decompression, bougie for sleeve sizing and OG tube for leak test. (Fig. 1A). Clinical risks of the current recommendation includes esophageal perforation, accidental stapling, and corkscrewing of the staple line.

Objective

The patient safety profile can be improved by reducing the above-mentioned steps. The aim of this study was to compare the current recommendation to a suction calibration system, ViSiGi 3D™, that performs all functions with one insertion, under a safe level of suction.

Materials/Methods

Unlike the regular bougie, anesthesiologist inserts ViSiGi 3D™ at the beginning of the case and only removes it once after stapling and leak testing is complete.



Comparison 1: Visual Confirmation of the Position


Primary sleeve gastrectomies were performed with a bougie or a ViSiGi 3D™ in a randomized, alternating order. Intraoperative pictures of the stomach prior to the first staple firing were captured in each case. (Fig. 3)

Comparison 2: Operating Time Savings

Times for completing stomach decompression, positioning, stapling, leak testing, and the total operating time, excluding complications, were obtained. (Fig. 4, Fig. 6)

Comparison 3: Staple Line Straightness

Three measurements were made at three locations on the inflated, excised gastric specimen: circumference, and distance from the greater curvature to both the anterior and posterior (Fig. 2). The variance of deviation from the midline at these locations was calculated for each specimen. (Fig. 5)



Comparison 4: Clinical Risk Reduction

The frequency of device movements for both groups were recorded. This record includes the number of intraluminal devices inserted as well as the number of times the surgeon inquired about the status/location of the devices. (Fig. 7)

Results

In a single-center randomized study, 26 patients were enrolled for a sleeve gastrectomy. The patient group consisted of 15 women and 11 men. The mean age was 36.8 years, ranging from 14 to 74 years. The overall BMI of patients averaged 45.3 kg/m².







Fig. 3 Side-by-side comparison of the stomach after device placement, immediately before the first stapling. In contrast to the bougie, ViSiGi 3D™ retained its position along the lesser curvature and provided a clearer delineation, hence reducing the risk of accidental tube stapling.



Conclusion

Comparison 1: The location of the bougie was difficult to ascertain, whereas ViSiGi 3D™ clearly delineated and indicated proper staple placement.

Comparison 2: ViSiGi 3D™ significantly decreases the time needed to complete each step of a sleeve gastrectomy, and therefore total operating time.

Comparison 3: ViSiGi 3D™ utilizes suction to maintain gastric placement, which equalizes tension on both sides of the stomach, thereby reducing corkscrewing.

Comparison 4: Fewer tube insertions may reduce esophageal damage and accidental tube stapling.

Acknowledgement and References

- R. J. Rosenthal. "International Sleeve Gastrectomy Expert Panel Consensus Statement: best practice guidelines based on experience of >12,000 cases". Surgery for Obesity and Related Diseases 8:1, 8-19
- Cost of treatment due to perforation, Oyasiji, tolutope. "Esophageal Perforation: Etiology, Outcome and Cost Analysis Over a Decade in a Community Teaching Hospital". SAGES 2012 Annual Meeting
- Abu-Gazala, Samir. "Nasogastric tube, temperature probe, and bougie stapling during bariatric surgery: a multicenter survey." Surgery for Obesity and Related Diseases. 8:5 (2011) 595-600
- Pictures in Fig. 5 are courtesy of St. Luke's Hospital - Allentown Campus
- This study was supported by Boehringer Laboratories, Inc.

Presented at:

19th World Congress of the International Federation for the Surgery of Obesity and metabolic disorders (IFSO2014).

5th International Consensus Summit for Sleeve Gastrectomy (ICSSG).

*August 26-30, 2014
Montreal, Canada*