



## The Infection Control Advantage!™

Boehringer is committed to improving healthcare by designing suction controls with features that aid infection control strategies. The Platinum Series Suction Regulators have distinct advantages to assist the effort to drive down the incidence of hospital-acquired infections.

### Self-Draining Gauge

The patented XL Gauge is uniquely designed to be exceptionally tolerant of flooding. In the event that overflow protection fails, the user can rely on the contaminants being evacuated effectively from the device.

Model 3800 Shown with optional integral trap bottle and antimicrobial cover



### Antimicrobial Knob Cover

Our optional silver-impregnated elasticized control knob cover is proven to provide a 10<sup>5</sup> reduction in bacteria in 24 hours. (MRSA test report available upon request.)

### Self-Clearing Technology™\*

Boehringer is able to offer a 12 year warranty on the 3800 Series Suction Regulators due in large part to the development of Self Clearing Technology™ that effectively protects the internal mechanism from contamination. Self Clearing Technology™ introduces airflow into the regulator body that sweeps aerosolized contaminants and fluids away from the internal mechanisms of the regulator. Keeping these elements clean and dry insures the proper operation of the device and reduces the risk of user contact with contaminated fluids.

\* Patent Pending

### UNLIMITED DISINFECTION AND STERILIZATION OPTIONS

All-metal materials of construction. Unlike plastic products, these devices are not vulnerable to cracking, discoloration and degradation. Compatible with Autoclave, EO, and any other process of choice.

**BOEHRINGER**®

Boehringer Laboratories, LLC  
500 East Washington Street  
P.O. Box 870  
Norristown, PA 19404

800-642-4945 [www.boehringerlabs.com](http://www.boehringerlabs.com)

Dwg 0000.031 Rev B

© 2007 Boehringer Laboratories, Inc. or its subsidiaries