

The Boehringer Model 7720 VCS was a specialty control made to facilitate Vacuum Assisted Venous Drainage. These devices were produced under a cooperative manufacturing arrangement between Boehringer Laboratories and Baxter Edwards in 1999 through 2000. These controls were originally built with a five (5) year warranty. Given the age of these units and the cooperative nature of their production, many of the parts needed to properly maintain them no longer exist. The lack of replacement parts was a major reason this unit was contemporized and the new model 3930 is offered as a direct replacement.



Obsolete Model 7720



Replacement 3930

Bearing in mind the critical nature of the procedures for which these instruments are used, and that their failure could have catastrophic results, facilities are strongly encouraged to consider updating to the Model 3930. Boehringer Laboratories will no longer accept Model 7720 units for Preventative Maintenance, Calibration, or Service and there are NO factory authorized Repair Centers for any of these devices.

The Model 3930 improved safety features over the original Model 7720 including - -

- Reduced 0-60mmHg working range in keeping with current clinical practices.
- Waste Anesthetic Scavenging in compliance with CDC/NIOSH mandates
- Independent Positive Pressure Relief Safety System
- Self-Clearing™ control systems to handle inadvertent fluid intrusion