

SUCTION REGULATOR SPECIFICATIONS

	Platinum Series 12 Year Warranty			Classic Series and Specialty Products 10 Year Warranty								
Model	3800	3804	3840	3700	3710	3714	3701	3704	3720	3740	3748	3750
Modes	OFF/ REG/ LINE	OFF/ CONTIN/ INTERMIT	OFF/ REG/ LINE	OFF/ REG/ LINE	OFF/ REG	OFF/ CONTIN/ INTERMIT	REG	OFF/ CONTIN/ INTERMIT	REG w/ exclusive Power Pulse feature	OFF/ REG/ LINE	OFF/ REG/ LINE	OFF/ REG
Range (mm Hg)	0-200 Full Line	0-200	0-760	0-200 Full Line	0-100	0-100	0-200	0-200	0-60	0-760	0-760	0-60 CM H ₂ O
Inlet/ Outlet	1/8" NPT											
Port Size	5 mm											
Spiking	+/- 10% at 200 mm Hg from full flow to zero											
Gauge Accuracy	+/- 5% FSO ANSI Class B											
Materials	Aluminum, Stainless Steel, Lexan, Buna Rubber, Silicone											
Dimensions (H x W)	5.75" x 2.75"	5.75" x 4.25"	5.75" x 2.75"	5.75" x 2.75"	5.75" x 2.75"	5.75" x 4.25"	5.75" x 2.5"	5.75" x 4.0"	5.75" x 3.0"	5.75" x 2.75"	5.75" x 2.75"	9.5" x 2.5"
Weight (lb)	1.3	1.65	1.3	1.3	1.3	1.65	1.2	1.65	1.4	1.3	1.3	1.85
Disinfection	Autoclave, Gas and Chemical Sterilants											

SUCTION REGULATOR SPECIFICATIONS

	Silver Series 5/3 Year Warranty			Classic Series 5/1 Year Warranty			
Model	7800	7804	7840	7700	7701	7702L	7740
Modes	OFF/REG/LINE	OFF/CONTIN/ INTERMIT	OFF/REG/LINE	OFF/REG/LINE	Regulate	OFF/CONTIN/ INTERMIT	OFF/REG/LINE
Range (mm Hg)	0-200 Full Line	0-200	0-760	0-200 Full Line	0-200	0-200	0-760
Inlet/ Outlet	1/8" NPT						
Port Size	5 mm						
Spiking	+/- 10% at 200 mm Hg from full flow to zero						
Accuracy	+/- 5% FSO ANSI Class B						
Materials	Aluminum, Stainless Steel, Lexan, Buna Rubber, Silicone						
Dimensions (H x W)	5.5" x 2.9"	5.5" x 4.25"	5.5" x 2.9"	5.5" x 2.9"	5.5" x 2.9"	5.5" x 4.25"	5.5" x 2.9"
Weight (lb)	1.4	1.8	1.4	1.4	1.3	2.6	1.4
Disinfection	Autoclave, Gas and Chemical Sterilants						