

# CareDry™ Regulator

*Maximizing the Performance of External Female Catheters*

## Standardized Care

- Pre-set suction ensures consistency of care and safety for all your patients.

## Fewer Leaks

- Up to 10x increased airflow accelerates drying time, helps prevent skin breakdown, and limits unnecessary barrel roll linen changes.

## Fewer Wounds

- A low, fixed pressure setting of 55mmHg provides a safe level of suction every time.

## Safe Access to Suction

- Does not require use of your existing suction regulators, thus maintaining suction availability for critical needs.



# The Advantages of Using the CareDry™ Regulator

Derived from a surgical suction device with more than a quarter-million uses, this device, along with the Power Take Off (PTO) ensures the results you need for your patients.



◀ The PTO maintains regulated suction availability for emergency use without capital purchases.

◀ With this configuration, the canister can be placed at any height with any length of tubing.

◀ The **CareDry™ Regulator** connects between the EFC and your cannister.

**BOEHRINGER**®

**PRODUCT ORDERING:**

CareDry™ Regulator	10 devices/ case	#7630
Power Take Off (PTO) Wall Adapter	1 ea.	#7622



Inservice Video

Caring for Lives through Innovation, Quality and Service

Boehringer Laboratories, LLC  
300 Thoms Drive Phoenixville, PA 19460  
800-642-4945 • [www.boehringerlabs.com](http://www.boehringerlabs.com)

0000.321 — CareDry™ Sales Sheet L321 Released Oct 2020