


CareDry® System

The Next Evolution in Urinary Management



The CareDry® System, with its novel antimicrobial sponge, is designed to free up nursing time while providing dryness and comfort for your patients.

- Up to 67% fewer product changes*
- 2x the absorbable surface area as leading competitors
- Does not rely on gravity for liquid collection



SCAN FOR DEMO



SCAN FOR WEBSITE


* indicated for up to 48 hrs of use, compared to 8-12 hr change interval of other female external catheters - data on file with Boehringer Laboratories

BOEHRINGER®

Your Best Solution for Urinary Management


24-hour Change Interval*

Antimicrobial sponge actively inhibits harmful bacteria formation on the sponge**, allowing for extended patient use.




100% Moisture-Wicking

Optimal pore size and foam structure with 5x the surface area of leading competitors provide rapid dryness and comfort for your patients.



Anatomically Conforming

Medical-grade foam protects against leaks in any position by contouring to the patient's unique body shape.



360-Degree Suction

Circumferential suction from a 4-channel core draws moisture away quickly from the patient and surrounding environment.

Product Ordering:

Description	Part Number	Ordering
CareDry® System	7660	20 units/box

* indicated for up to 48 hrs of use, compared to 8-12 hr change interval of other female external catheters

** data on file with Boehringer Laboratories